


**Life Fitness Models 9100, 9500HR, 9700HR, and 9700HR w/Decline Treadmills**  
**DIAGNOSTICS: SYSTEM TEST – TELEMETRY MODE**

---



Upon initial entry to this test, a message concerning the telemetry system on/off will occur. This test will allow the user to test the telemetry heart rate system. If the telemetry heart rate system is OFF, no heart rate value will be given.

If the system is enabled, a telemetry heart rate signal can be provided and a heart rate value will be shown. The signal will also generate the 'ENTER' led to flash at the heart rate pulse frequency.

Press the 'CLEAR' key to exit the Telemetry Heart Rate test and return to the System Test Menu.